

MICRO SURGERY

instruments





Street # 7, ShahabPura Road, Sialkot 51310 - Pakistan

Ph: +92-52-3553747 Fax: +92-52-3552417

info@lebord.com www.lebord.com



Paton Needle Holder

LBO-101



Heavy11mm serrated jaws with lengthwise groove. Flat serrated handle Overall length 119mm, 4.7 inches.



Catalano Needle Holder-Curved with Lock

LBO-102



Fine 9mm moderately curved jaws with flattened bottom. Large round knurled handle with lock. Overall length 135mm, 5.3 inches.



Castroviejo Needle Holder-Curved

LBO-103 LBO-104



Medium 11mm gently curved jaws. Flat serrated handle. Overall length 137mm, 5.4 inches.

LBO-103 with Lock
LBO-104 without Lock





Castroviejo Needle Holder-Straight with Lock

LBO-105



Azar Needle Holder-Curved with Lock

LBO-106



Castroviejo Needle Holder-Straight with Lock

LBO-107



Castroviejo Needle Holder-Straight with Lock

LBO-108





Castroviejo Needle Holder-Straight with Lock

LBO-109



Castroviejo Needle Holder-Curved with Lock

LBO-110



O' Gawa Needle Holder-without Lock Straight

LBO-111



Delicate 8mm jaws. with well rounded edges to prevent cutting suture. Overall length 130mm, 5.1 inches.



Needle Holder-Curved

LBO-112 LBO-113



Delicate 6mm moderately curved jaws. Small round knurled handle. Overall length 128mm, 5.0 inches.

LBO-112 Without Lock
LBO-113 With Lock





Troutman Needle Holder-Curved

LBO-114 LBO-115



Medium 10mm gently curved jaws. Medium round "horizontally" serrated handle. Overall length 122mm, 4.8 inches.

LBO-114	Without Lock
LBO-115	With Lock



Sinskey Needle Holder-Straight with Lock

LBO-116



Delicate 7mm jaws. Smooth flat handle with lock. Overall length 111mm, 4.4 inches.



Castroviejo Needle Holder

LBO-117



Straight 9mm medium jaws. Round knurled handle with lock. Overall length 134mm, 5.3 inches.





Kalt Needle Holder





Heavy 12mm straight serrated jaws with tungsten carbide inserts. Flat serrated locking handle. Overall length 141mm, 5.6 inches.

Kalt Needle Holder





Heavy 12mm straight jaws with edges of top jaw serrated. Flat serrated locking handle . Overall length 141mm, 5.6 inches.

Kalt Needle Holder-Original Size Jaw

LBO-120



 $\label{lem:medium 10mm} \mbox{Medium 10mm straight serrated jaws. Flat serrated Locking handle} \mbox{Overall length 141mm, 5.6 inches.}$

Kalt Needle Holder





Medium 10mm straight serrated jaws with jaws 20% finer than original model. Flat serrated Locking handle . Overall length 137mm, 5.4 inches.





Sharp pointed tips. 16mm from mid-screw to tip. Wide serrated handle Overall length 132mm, 5.2 inches.

Westcott Curved Tenotomy Scissors



Rounded Westcott blades with blunt tips. 19mm from mid-screw to tip. Flat serrated handle. Overall length 115mm, 4.5 inches.

LBO-123 Right LBO-124 Left

Shepard-Westcott Curved Tenotomy Scissors LBO-125

Right model. Westcott blades with blunt tips. Upper blade serrated. 21mm from mid-screw to tip. Flat serrated handle . Overall length 123mm, 4.8 inches.

Fine Straight Stitch Scissors LBO-126

Sharp pointed tips. 17mm from mid-screw to tip. Flat serrated handle Overall length 112mm, 4.4 inches.





LBO-127

Castroviejo

Universal Corneal Scissors curved blunt tips, large blades



LBO-128

Castroviejo

Universal Corneal Scissors strongly curved blunt tips, large blades



LBO-129

Castroviejo

Keratoplasty Scissors angled to side blunt tips, medium blades



LBO-130

Castroviejo

Conjunctival Scissors curved blunt tips, small blades



LBO-131 DeWecker Iris Scissors



Blunt and sharp points. 12mm from bend to tip. Flat serrated squeeze action handle. Overall length 114mm, 2.1 inches.



Miniature Corneal Scissors




Delicate corneal scissors for use during corneal transplants to complete the corneal trephination. 9mm blades with lower blade extended 0.5mm. Overall length 100mm, 3.9 inches.

Katzin Curved Corneal Transplant Scissors



Miniature model with strongly curved blunt blades of uneven lengths. lower blade 12mm, upper blade 11mm (from mid-screw to tip). Flat serrated handle. Overall length 108mm, 4.2 inches.





Gills-Welsh Vannas Angled Micro Scissors

LBO-138



11mm blades with sharp pointed tips. 13mm from mid-screw to tip. Serrated handle Overall length 86mm, 3.4 inches.

Vannas Angled Capsulotomy Scissors

LBO-139



6mm blades with sharp pointed tips. 9mm from mid-screw to tip. Serrated handle. Overall length 81mm, 3.3 inches.

Capsule Scissors 10mm Angled Blades

LBO-140



Delicate 10mm blades angled \$45 for performing capsulotomies and snipping capsule tags, iris, and other Structures in the anterior chamber. Overall length 110mm, 4.3 inches.

Capsule Scissors 10mm Curved Blades

LBO-141



Delicate 10mm blades angled 45. Blades gently curved to offer exception1 control in the anterior chamber when performing anterior capsulotomies and snipping capsule tags. Overall length 110mm, 4.3 inches.



Vannas Straight Scissors Vannas Curved Scissors

LBO-142





LBO-143



Sharp pointed tips. 7mm from mid-screw to tip. Smooth handle . Overall length 83mm, 3.2 inches.

Vannas Straight Capsulotomy Scissors Vannas Curved Capsulotomy Scissors

LBO-144





LBO-145



6 mm blades with sharp pointed tips. 8 mm from mid-screw to tip. Serrated Overall length $84 mm,\,3.3$ inches.

Vannas Angled Capsulotomy Scissors

LBO-146



6mm blades with sharp pointed tips. 9mm from mid-screw to tip. Serrated handle . Overall length 81mm, 3.1 inches.

Gills-Welsh Vannas Angled Micro Scissors

LBO-147




5mm blades with sharp pointed tips. 8mm from mid-screw to tip. Serrated handle . Overall length 80mm, 3.1 inches.



Wilmer Angled Conjunctival and Utility Scissors







Elsching Fixation Forceps

LBO-152





LBO-153

Straight shafts with 2.0mm, 1 x 2 Elschnig style teeth. Serrated handle . Overall length 109mm, 4.3 inches.

LBO-152 Narrow Wide

Wide handle

Oshre Conjunctival Fixation Forceps

LBO-154





Straight shafts with 5.5mm tying platforms and 45 0.3mm teeth. Wide serrated handle Overall length 108mm, 4.3 inches.

Bracken-Farkas Corneal Forceps with Neuman Modification

LBO-155





Straight shafts with delicate 0.12mm, 1 \times 2 teeth. Stabilization pin prevents twisting of tissue. Serrated handle . Overall length 107mm, 4.2 inches.

Iowa State Fixation Forceps

LBO-156





Straight shafts with 1 x 2 teeth. Flat handle . Overall length 86mm, 3.4 inches.



Osher 0.12mm Forceps - 45 Angled Teeth

LBO-157





The teeth of this forceps permit precise control of the wound lips. 3mm tying platform handles the finest suture material. 0.12mm, 1×2 teeth. Overall length 104mm, 4.1 inches.

Osher 0.12mm Forceps-Elschnig Teeth

LBO-158





The Elschnig teeth of this forceps are ideal for grasping sclera during suturing. Delicate tying platform handle the finest suture material. 0.12mm, 1×2 teeth. Overall length 104mm, 4.1 inches.

Green Fixation Forceps

LBO-159





Narrow 12mm jaws with fine teeth. Serrated handle Overall length 100mm, 3.9 inches.

Kremer Fixation Forceps

LBO-160





Kremer Angled Fixation Forceps

LBO-161





Two point fixation, tips separated by 13mm. Each tip has 0.12mm, 1 x 2 teeth. Wide serrated handle with thumb catch lock.

Overall length 105mm, 4.1 inches.

LBO-160 Straight



Snellen Entropion Forceps-Small

LBO-162



Solid lower plat with outside width 28mm. Open upper plat. Serrated handle wit locking thumb screw. Overall length 88mm, 3.5 inches.

LBO-163



LBO-162 Left LBO-163 Right

Wies Chalazion Forcep

LBO-164



Solid lower plate with outside width 15mm. Open serrated upper plate with inside dimensions 9mm x 12mm. Flattened top permits placement deep into the fornix. Serrated handle with locking thumb screw. Overall length 94mm, 3.7 inches.

Demarres Chalazion Forceps

LBO-165



Overall solid 16mm wide lower plate. Open upper plate with inside dimensions $12mm \times 14mm$. Cross action serrated handle . Overall length 96mm, 3.8 inches.

Demarres Chalazion Forceps

LBO-166



Oval solid 24mm wide lower plate. Open upper plate with inside dimensions 12mm x 20mm. Serrated handle with locking thumb screw . Overall length 92mm, 3.6 inches.



Francis Chalazion Forceps

LBO-167



Oval solid 17mm wide lower plate. Open upper plate with dimension 12mm x 14mm. Serrated handle with locking thumb screw . Overall length 97mm, 3.8 inches.

Heath Chalazion Forceps

LBO-168



Oval solid 16mm wide lower plate. Open upper plate with dimension 12mm x 14mm. Serrated handle with locking thumb screw . Overall length 96mm, 3.8 inches.

Hunt Chalazion Forceps

LBO-169



Round solid 15mm wide lower plate. Open upper plate with inside diameter 12mm. Serrated handle with locking thumb screw . Overall length 97mm, 3.8 inches.

Lambert Chalazion Forceps-Small

LBO-170



Round solid 10mm wide lower plate. Open upper plate with inside diameter 8mm. Serrated handle with locking thumb screw . Overall length 92mm, 3.6 inches.

Baird Chalazion Forceps

LBO-171



Oval solid 15mm wide lower plate. Open upper plate with inside dimensions 11mm x 8mm. Serrated handle with locking thumb screw . Overall length 97mm, 3.8 inches.



Osher Capsule-Foreign Body Forceps

LBO-172





Fine 3.0mm platforms Angled 60 10mm from angle to tip. Serrated handle. Overall length 80mm, 3.1 inches.

Gills-Welsh Capsule Forceps, Colibri Style

LBO-173





Smooth Colibri tip used to grasp and subsequently retract the capsule during nucleus delivery. Serrated handle . Overall length 73mm, 2.9 inches.

Harms Corneal Utility Forceps, Colibri Style

LBO-174





0.12mm, 1 x 2 teeth set at 45 with 4.5mm tying platform. Serrated handle Overall length 70mm, 2.7 inches.

Cohan Corneal Utility Forceps, Colibri Style

LBO-175





0.12mm, 1 x 2 teeth set at 45 with 4.5mm tying platform. Serrated handle. Overall length 60mm, 2.4 inches.

Capsulorhexis Forceps

LBO-176





The ultra-fine tips are designed with an iris stop and angled platform tips for grasping delicate capsule. Overall length 85mm, 3.3 inches.

Ultra Capsulorhexis Forceps

LBO-177





Fine angled tips ideal for grasping delicate capsule. Angled 12mm shafts with iris stop allows ease of movement in the anterior chamber. Flat handle. Overall length 80mm, 3.2 inches.





Pierse Forceps-Colibri

LBO-178





Pierse tips with 6mm platforms and 2.5mm angled tips. Flat serrated handle. Overall length 86mm, 3.3 inches.

Pierse Forceps-Straight

LBO-179





Straight shafts with 6mm platforms and Pierse tips. Flat serrated handle. Overall length 85mm, 3.3 inches.

Moody Curved Fixation Forceps

LBO-180





Curved shafts with 0.5mm, 1 x 2 teeth. Serrated handle with thumb catch lock . Overall length 108mm, 4.3 inches.

Great Fixation Forceps

LBO-181









THUMB CATCH LOCK

LBO-181	Narrow, 3.5mm jaws	Without
LBO-182	Wide, 4.5mm jaws	Without
LBO-183	Wide, 4.5mm jaws	With

Jaws with teeth. Serrated handle . Overall length 108mm, 4.3 inches.



Maumenee Corneal Forceps, Colibri style

LBO-184





LBO-185



5mm tying platform, serrated cross action handle Overall length 121mm, 4.8 inches.

	TEETH
LBO-184	0.10mm, 1 x 2 set at 45
LBO-185	0.10mm, 1 x 2 set at 90

Corneal Forceps, Colibri Style (Double)

LBO-186





Two 2.75mm long tips separated 1mm with 0.4mm Pierse tips. Serrated handle. Overall length 72mm, 2.8 inches.

Polack Double Corneal Forceps, Colibri Style

LBO-187





Two 2.75mm long tips separated 1mm with 1 x 2 teeth. Serrated handle. Overall length 71mm, 2.8 inches.

Lieberman Suturing Forceps

LBO-188





Colibri style tips. Double toothed tips with 0.10mm teeth spaced 1mm apart and 7mm tying platforms. Wide serrated handle . Overall length 109mm, 4.3 inches.



Cilia Forceps

LBO-189





Oval solid 15mm wide lower plate. Open upper plate with inside dimensions $11mm \times 8mm$. Serrated handle with locking thumb screw . Overall length 97mm, 3.8 inches.

Cilia Forceps

LBO-190





Narrow smooth jaws. Serrated handle Overall length 83mm, 3.3 inches.

Barraquer Cilia Forceps

LBO-191





Also referred to as McLean-Putney suture forceps. Slightly curved 6mm long platforms. Serrated handle . Overall length 109mm, 4.3 inches.

Ziegler Cilia Forceps

LBO-192





Slender 4.5mm long rounded platforms. Serrated handle Overall length 91mm, 3.6 inches.

Berke Ptosis Forceps

LBO-193



longitudinal serrations. Slide lock and serrated handle.



O'Connor-Elsching Fixation Forceps

LBO-194





Straight shafts with heavy Elschnig teeth. Flat serrated handle with slide style lock . Overall length 109 mm, 4.3 inches.

Schaaf Foreign Body Forceps

LBO-195





Fine pointed grooved jaws with smooth curve from tip to finger grip style handle. Overall length 105mm, 4.1 inches.

Rappazzo Intraocular Foreign Body Forceps

LBO-196



LBO-197



LBO-198



LBO-296 2.0mm smooth
LBO-297 1.5mm Electro Dusted

LBO-298 1.5mm Serrated

Inner jaw space 2.0mm-45 tip angulation. 20 gauge (.89mm) diameter shaft and stationary distal jaw. Lightweight squeeze handle. Overall length 163mm, 6.4 inches.


Rizzuti Superior Rectus Forceps

LBO-199





Shaft angled 45, 7mm from tip to bend. Tips curved with 1 x 2 teeth. Serrated cross action handle . Overall length 97mm, 3.8 inches.








Curved Dressing Forceps Serrations





Gently curved tips. Serrated handle . Overall length


LBO-203



		OVERALL LENGTH
LBO-202	Delicate	99mm,3.9 inches
LBO-203	heavy	98mm,3.9 inches

Curved Dressing Forceps Serrations

LBO-204





Delicate strong curved serrated tips. Serrated handle . Overall lenght 97 mm, 3.8 inches.

Bishop-Harmon Dressing Forceps

LBO-205






Cross serrated tips.Flat handle .Overall lenght 100 mm, 3.9 inches.

Sinsky-Mcpherson Angled Tying Forceps

LBO-206





Angled 45 Shafts 7.5mm From Tip To Bend With 4mm Tying Platforms And Irrigator. Serrated Handle. Uses 40cm Of Silicone Tubing.



Knolle-Shepard Lens Forceps

LBO-207



Curved shafts with 2.0mm wide serrated tips and 0.4mm pin in upper jaw. Serrated cross action handle. Overall length 121mm, 4.8 inches.



Cozean-McPherson Angled Lens Forceps

LBO-208





Angled 45 shafts 7mm from tip to bend with grooved 2mm platform for grasping lens haptic. Serrated McPherson style handle. Overall length 80mm, 3.1 inches.

Cozean Angled Lens Forceps

Lehner IOL Forceps

LBO-209





12 mm jaws with 8 mm angled tips. Serrated cross action handle. Overall length $98 mm,\, 3.9$ inches.

LBO-210





Gently curved serrated tips make this lens forceps ideal for tunnel incisions. Round handle. Shaft angled 11mm from tip. Overall length 103mm, 4.1 inches.

Clayman Curved Tying Forceps

LBO-211





Curved shafts and tips with 6mm tying platforms. Wide smooth handle. Overall length 100mm, 3.9 inches.



LBO-212



Nordan Colibri Forceps



0.12mm 1 x 2 teeth with 5mm platforms and 2.0mm angled tips. Round knurled handle. Overall length 113mm, 4.5 inches.

Kirby Corneoscleral Forceps.

LBO-213





 $1\ x\ 2$ teeth, cross serrated 5mm 60 angled platforms. Round knurled handle . Overall length 107mm, 4.2 inches.

Catalano Corneal Scieral Forceps

LBO-214





1 x 2 teeth and 5.5mm tying platforms. Round knurled handle. Overall length 113mm, 4.4 inches.

TEETH LBO-214 0.12mm **LBO-215** 0.3mm

Nordan Colibri Forceps

LBO-216





Stainless steel model of the popular Nordan titanium Colibri forceps. tying platform with 0.12mm teeth for grasping delicate structures. Flattened round handle Overall length 104mm, 4.1 inches.

Katzin-Barraquer Corneal Utility Forceps, Colibri Style

LBO-217






0.4mm, 1×2 teeth and 7mm tying platform. Serrated handle Overall length 73mm, 2.9 inches.





LBO-218





LBO-219

Straight shafts with 1mm, 45 angled tips. Fine serrations. Smooth jeweler's forceps style handle. Overall length 109mm, 4.3 inches.

Hydro Insertion Forceps

LBO-220



Forceps designed for insertion of foldable intraocular lenses. Unique cross-action design allows for easy insertion of most foldable lens style. Overall length 119mm, 4.7 inches.

LBO-221





This instruments is ideal for cracking the lens nucleus as well as insertion of foldable intraocular lenses. Overall length 90mm, 3.5 inches.

Beaupre Cilia Forceps

LBO-222






1mm smooth 45 angled jaws. Cross action serrated handle. Overall length 117mm, 4.6 inches.



LBO-250 to LBO-253





Knapp

Lacrimal Retractor 4 prongs

A: blunt

B: sharp

LBO-259

LBO-258

Helveston

Tissue Retractor thin curved blade

- A: "Barbie" 7mm wide
- **B:** "Big Barbie" 9mm wide
- C: "Great Big Barbie" 11mm wide concave blade







20mm x 35mm & 24mm x 35mm curved ends Clear Lucite Overall Length 93mm 3.7 inches



Curved over length & width 20mm & 23mm wide blades. Overall length 109mm, 4.3 inches.



Thin ribbon retractor with a gentle "s" curve. Double ended 16mm wide & 25mm wide ends. Overall length 182mm, 7.2 inches.



Thin ribbon retractor with a gentle "s" curve. Double ended 10mm wide & 14mm wide ends. Overall length 180mm, 7.1 inches.



Smooth 14mm by 56mm blades with 4.5 wide notch in blade. Flat serrated handle. Overall length 148mm, 5.8 inches.



LBO-265

Stevenson

Lacrimal Sac Retractor 3 x 3 blunt prongs. Adjustable.



LBO-266

Stevenson

Lacrimal Sac Retractor solid blades with serrated edge adjustable



Kennerdell Spatulas



Set of 10 blades (two each of five sizes). Malleable Construction allows easy bending.

LBO-267 to LBO-271




	SIZE (mm)	OVERALL LENGTH
LBO-267	6 x 49	92mm, 3.6 inches
LBO-268	6 x 74	117mm, 4.6 inches
LBO-269	14 x 62	105mm, 4.1 inches
LBO-270	18 x 60	104mm, 4.1 inches
LBO-271	20 x 66	110mm, 4.3 inches





Barraquer Eye Speculum

Wire 14mm blades. Blade spread 18mm. Overall length 38mm, 1.5 inches.




Kratz-Barraquer Eye Speculum

Open 14mm blades. Blade spread 18mm. Overall length 41mm, 1.6 inches.



Temporal Approach Lid Speculum-Wire Blade

Angled for temporal approach. Wire 14mm rounded blades for patient comfort. Blade spread 18mm. Overall length 38mm, 1.5 inches.




Temporal Approach Lid Speculum-Open Blade

Angled for temporal approach. Open 14mm rounded blades for patient comfort. Blade spread 18mm. Overall length 41mm, 1.6 inches.



Barraquer Eye Speculum

Open 14mm blades of heavy 1mm wire. Blade spread 20mm. Overall length 40mm, 1.6 inches.



Kratz-Nichamin Speculum

Open wire style 16mm blade with a wing tip twist giving the surgeon better on-axis accessibility. Overall length 45mm, 1.8 inches.





Barraquer Eye Speculum

Solid 14mm blades. Blade spread 18mm. Overall length 41mm, 1.6 inches.




Barraquer Eye Speculum-Pediatric

Open 11mm Blade. Blade spread 19mm Overall length 34mm, 1.3 inshes.



Temporal Approach Lid Speculum-Solid Blade

Angled for temporal approach. Solid 14mm blades. Blade spread 18mm. Overall length 41mm, 1.6 inches




Kartz-Barraquer Eye Speculum-Pediatric

Open 10mm blades. Blade spread 19mm. Overall length 35mm, 1.4 inches.



Barrapuer Pediatric Speculum-Solid Blade


Solid 6mm solid blade speculum for newborn or pediatric use. Overall length 30mm, 1.2 inches



Barraquer Pediatric Speculum-Wire Blade

Closed wire 5mm blade with a wing tip twist giving th surgeon better on-axis accessibility. Overall length 30mm, 1.2 inches.






Floyd-Barraquer Eye Speculum

Open 14mm blades with 6mm lifting wings. Blade spread 18mm. Overall length 40mm, 1.6 inches.



Sauer Eye Speculum-Infant

Solid 8mm blades. Blade spread 20mm. Overall length 35mm, 1.4 inches.



Simcoe Eye Speculum

Open 14mm blades with 10mm lifting wings. Lightweight wire. Blade spread 18mm. Overall length 40mm, 1.6 inches.



Saunders Eye Speculum-Premature Infant

Solid 4mm blades. Locking thumb screw. Blade spread 14mm. Holes in bar decrease weight of speculum.

Overall length 35mm, 1.4 inches.



Pierse Eye Speculum

Open 11mm blades. Blade spread 20mm. Overall length 30mm 1.2 inches.




Alfonso Eye Speculum-Newborn

open 5mm blades. Blade spread 27mm. Overall length 45mm, 1.8 inches.
















Aspirating Speculum with open blades
For unobstructed access to the cataract incision with phaco or I/A tip.
Supplied with silicone tubing and luer-lock adaptor.









Mellinger Eye Speculum

Fenestrated 16mm blades. Self-locking Mechanism. Blade spread 27mm. Overall length 70mm, 2.8 inches.



16mm open wire blades. Self locking mechanism. Blade spread 27mm. Overall length 70mm, 2.8 inches.



MODEL	BLADES	BLADES SPREAD	OVERALL LENGTH
R58-4908 Infant	7mm	24mm	52mm, 2.0 inches
R58-4909 Pediatric	10mm	22mm	55mm, 2.5 inches



Fenestrated 14mm blades. Self-locking mechanism. Blade spread 24mm. Overall length 57mm, 2.2 inches.



Mellinger Eye Speculum

Solid 16mm blades. Self-locking mechanism. Blade spread 27mm. Overall length 70mm, 2.8 inches.



Fenestrated 14mm blades. Wide extension of blade for protection of lids. Self-locking mechanism. Blade spread 23mm. Overall length 55mm, 2.2 inches.







LBO-304

Keizer-Lancaster Eye Speculum

Solid 15mm blades with 8mm narrowed portion. Screw-controlled locking mechanism with stabilizing disk. Blade spread 30mm. Overall length 82mm, 3.2 inches.



Open 13mm blades. Screw-controlled locking mechanism. Blade spread 60mm. Overall length 85mm, 3.3 inches.



Lester-Burch Eye Speculum

Solid 21mm blades. Screw-controlled locking mechanism. Blade spread 32mm. Overall length 80mm, 3.1 inches.



Maumenee-Park Speculum

Fenestrated 14mm blades and canthus bar. Adjustable blade locking screws on outside of branches. Blade spread 36mm. Overall length 80mm, 3.1 inches.




Osher Lid Spculum

Open wire style blades allow complete manipulation of the phaco tip inside the eye. Lid control is excellent, exposure is maximized, and the likelihood of positive pressure is reduced. Overall length 150mm, 5.9 inches.





Fenestrated 15mm blades and canthus bar. Adjustable blade locking screws on outside of branches. Blade spread 36mm. Overall length 82mm, 3.2 inches.



Spring-controlled spread of blades with locking mechanism.

BLADES	BLADES SPREAD	OVERALL LENGTH
14mm	52mm	86mm, 3.4 inches



Solid 14mm blades and canthus bar. Adjustable blade locking screws on outside of branches. Blade spread 36mm. Overall length 82mm, 3.2 inches.



Williams Eye Speculum

Spring-controlled spread of blades with locking mechanism.

BLADES	BLADES SPREAD	OVERALL LENGTH
11mm	41mm	69mm, 2.7 inches

LBO-313



McPherson Speculum

Solid 14mm free-swinging blades, suture posts and stabilizing disk. Blade spread 30mm. Overall length 82mm, 3.2 inches.



Solid 15mm blades. Screw controlled locking mechanism with stabilizing disk. Blade spread 35mm. Overall length 82mm, 3.2 inches.



Cannulas

Simcoe Anterior Chamber Receiving Cannulas





25 gauge with 35 curved shaft and beveled tip. Facilities exit of retaining wire. Overall length 15.5mm, 0.6 inches.

McIntvre Straight Lacrimal Cannula



LBO-327

23 gauge cannula with blunt tip and side port. Straight Shaft Overall length 46mm, 1.8 inches.

Simcoe Suturing Cannulas






30 gauge cannula with beveled sharp tip. Hole for intraocular suturing. Overall length 17mm, 0.7 inches.

McIntyre Curved Lacrimal Cannula





23 gauge cannula with blunt tip and side port. Overall length 46mm, 1.8 inches.

Lambert 33 Gauge Subretinal Cannulas



LBO-329






LBO-330

Designed for aspiration of fluid from subretinal space. 33 gauge tip. Overall lenght 35mm, 1.4 inches.

LBO-329 blunt sharp **LBO-330**

Lacrimal Cannula



LBO-331

23 gauge cannula with blunt tip, straight shaft, Overall length 25mm, 1.0 inches.

Atkinson Retrobulbar Cannulas



Beveled sharp rounded point, Overall length 35mm, 1.4 inches.

LBO-332 23 Gauge

LBO-333 25 Gauge

R59-4933 Shahinian Lacrimal Cannula





23 gauge cannula with blunt tapered tip and side port. Straight shaft . Overall length 11mm, 0.4 inches.

Cannulas



Lacrimal Cannula LBO-335





23 gauge cannula with reinforced shaft, straight 10mm tip, Overall length 31mm, 1.2 inches.

LBO-337 Simcoe Posterior Capsule Polisher





23 gauge cannula with sandblasted tip. Shaft angled45, 14mm from tip to bend. Overall length 16mm, 0.6 inches.

Lacrimal Cannula LBO-336





23gauge cannula with reinforced shaft, curved 10mm tip,
Overall length 31mm, 1.2 inches.

LBO-338 Thornton Aspirating Capsule Polisher





Aspirating 27 gauge cannula with "U" shape tip permits aspiration of 12 o'clock cortex. The tip is sandblasted to act as a capsule polisher. Overall length 23mm, 0.9 inches.

Weil Lacrimal Cannula LBO-339





21gauge cannula with blunt tip and top port. Curved shaft . Overall length 34mm, 1.3 inches.

Gills LBO-340





Aspirating Cannula smooth blunt tip end opening angled, 25 gauge

Tip Size

A: 5mm tip

B: 7mm tip

C: 8mm tip

D: 10mm tip

Drews Capsule Polisher LBO-341





25 gauge cannula. 2.5mm diameter angled ring with port on inside at 6 o'clock. Sandblasted posterior surface, smooth anterior surface. Overall length 23mm, 0.9 inches.






LBO-342


Gimbe

Irrigating Cannula blunt tip, angled to side end opening 26 gauge



Cannulas





A: right
B: straight

C: left



LBO-347



Kelman

Cystotome sharp point and sharp cutting edge 22 gauge

LBO-348






Kratz

Cystotome sharp point and sharp cutting edge formed, 22 gauge

LBO-344 Simcoe Double-Barreled Irrigating-Aspirating Unit LBO-345 Simcoe Reverse Double-Barreled Irrigating-Aspirating Unit



LBO-344 Irrigation

Irrigation through main hub, aspiration through silicone tubing.

LBO-345

Irrigation through silicone tubing, aspiration through main hub



LBO-350

Silicone Tubing

For I/A Cannulas tubing and adaptor (Tubing only,12")




Scobee Oblique Muscle Hook

LBO-351



Curved shaft with flat 4.5mm hook. Flat serrated handle . Overall length 134mm, 5.3 inches.

LBO-352



7 mm flattened hook. Flat serrated handle . Overall length 125 mm, 4.9 inches.

Culler Rectus Muscle Hook

LBO-353



10 mm flattened hook with angled tip. Flat serrated handle . Overall length 122mm, 4.8 inches.


Jamson Muscle Hook

LBO-354



2.0 mm bulbous tip with 6mm flat portion. Flat serrated handle . Overall length 136mm, 5.4 inches.

LBO-355



Small bulbous tip with flat sides. Flat serrated handle . Overall length 132mm, 5.1 inches.



Green Strabismus Hook

LBO-356



Flat 7mm wide hook. Flat serrated handle. Overall length 130mm, 5.1 inches.



Graefe Strabismus Hook

LBO-357



Flattened hook. Flat serrated handle Choose from three sizes.

HOOK OVERALL LENGTH

A: # 1 size 8mm 135mm, 5.3 inches

B: # 2 size 10mm 140mm, 5.5 inches

C: #3 size 10mm 143mm, 5.7 inches



Greafe Strabismus Hook-Non-Magnetic

LBO-358



#2 Size 10mm hook. Flat serrated handle. Overall length 139mm, 5.5 inches.



LBO-359



3mm forked angled blades. Blade spread 7mm. Flat serrated handle Overall length 123mm, 4.8 inches.



Stevens Curved Tenotomy Hook

LBO-360



6 mm curved hook. Flat serrated handle . Overall length 119 mm, 4.7 inches.



Maumenee Iris Hook

LBO-361



Sharp, 1.5mm hook. Flat serrated handle . Overall length 123mm, 4.8 inches.

Shepard Iris Hook

LBO-362



Curved 1.0mm hook. Flat serrated handle . Overall length 123mm, 4.8 inches.

Bonn Iris Hook

LBO-363



Curved 0.3mm hook. Flat serrated handle. Overall length 123mm, 4.8 inches.

Hirschman Iris Hook

LBO-364



"H" shaped 0.75mm hook angled 45, 9mm from tip to bend. Round knurled handle. Overall length 122mm, 4.8 inches.

Kanses Nucleotome Bisection Spatula

LBO-365



This semi-sharp blade is designed for use in the left hand in conjunction with the Kansas Nucleus Support Spatula to split the lens nucleus as described by the "Phacofracture" technique. Round knurled handle. Overall length 115mm, 4.5 inches.

Hardesty Tenotomy Hook

LBO-366



5mm curved hook with marking at 14mm and 20mm from tip. Flat serrated handle . Overall length 120mm, 4.7 inches.

Smith Expressor Hook

LBO-367



11.5mm curved hook with marking at 14mm and 20mm from tip. Flat serrated handle . Overall length 150mm, 5.9 inches.

Gass Retinal Detachment Hook

LBO-368



Flattened 13mm "Graefe type" hook, with 1.5mm oval hole. Flat serrated handle . Overall length 143mm, 5.6 inches.

Lewicky Hook

LBO-369



0.2mm, 90° angled tip with curved 9mm shaft. Flat smooth handle . Overall length 115mm, 4.5 inches.

Azar Lens Manipulation Hook

LBO-370



0.2mm, 90 angled tip. Flat smooth handle . Overall length 120mm, 4.7 inches.



Wilder Lens Loop





4mm wide x 15mm long loop with serrations on top of loop. Flat serrated handle Overall length 135mm, 5.3 inches.

Lewis Lens Loop

LBO-372



6mm wide x 7mm long loop. Flat serrated handle . Overall length 135mm, 5.3 inches.

Gills-Welsh Modified Lens Loop

LBO-373



6mm wide x 7mm long loop. Flat serrated handle . Overall length 135mm, 5.3 inches.

Gills-Welsh Lens Loop with Olive Tip

LBO-374



5mm wide x 7mm long loop ultra thin loop. Olive tip on opposite end. Flat serrated handle . Overall length 121mm, 4.8 inches.





Gills-Welsh Morrison Lens Loop with Olive Tip

LBO-375




5mm wide x 7mm long loop ultra thin loop with 20 additional curvature. Olive tip on opposite end. Flat serrated handle . Overall length 121mm, 4.8 inches.











Kallan Nucleus Rotator

LBO-382



Multi-purpose manipulator with broad tip provides stability and control of the nucleus during phacoemulsification. The instruments can also be introduced through the side port incision during implantation of the IOL guide the lens down into the capsular bag. Round knurled handle. Overall length 125mm, 5.0 inches.

MINI O'Gawa Mini Lens Dialer

LBO-383



Nucleus manipulator facilitates lens rotation during phacoemulsification. Overall length 120mm, 4.7 inches.

Nichamin Double Ended Nucleus Chopper-Spatula

LBO-384



Wedge-shaped chopper designed to allow nucleus chopping in three directions.

Rounded tip highly polished to protect the posterior capsule. 0.5mm spatula end is used to manipulate softer nuclei. Round knurled handle. Overall length 140mm, 4.6 inches.



Minami "M" Hook

LBO-385



This tear dorp-shaped tip with multi-angled shaft is ideal for the "Miracle Multiple-Cut" technique and is used for dividing and manipulating nucleus material from all directions in the capsule. Round knurled handle.

Overall length 117mm, 4.6 inches.

Nagahara Nucleus Chopper




Nucleus divider with 1.25mm tip and blade offset at a 60 angle is ideal for performing the "Phaco Chop" technique. Round knurled handle. Overall length 114mm, 4.5 inches.

A: Left B: Right

LBO-386









Thornton Phacoemulsification Spatula

LBO-392



Short handled spatula for use through the side port stab incision. The malleable silver 0.5mm tip can be shaped to the surgeon's preference for maximum control. Round serrated handle. Overall length 86mm, 3.4 inches.



LBO-393



Double ended with needle point and 5mm by 14mm thin curved blade. Flat serrated handle . Overall length 138mm, 5.4 inches.



LBO-394



Flared 3mm wide by 14mm long platinum blade on 35mm shaft. Flat serrated handle . Overall length 140mm, 5.5 inches.



LBO-395



mmmmm

LBO-396



Drysdale Nucleus Manipulator

LBO-397



Flattened paddle tip with generous surface area is ideal for manipulating tissue, rotating, and cracking the nucleus. Round knurled handle. Overall length 110mm, 4.3 inches.

Olk Angled Retinal Spatula

LBO-398



Triangular 0.5mm angled tip with 4mm flat surface. Round knurled handle with dull finish Overall length 125mm, 4.9 inches. Tip Angles 115

Barrett Micro Mushroom Manipulator

LBO-399



Ideal for nucleus manipulator during phacoemulsification. Also used for manipulating and rotating IOLs within the anterior chamber. Round knurled handle. Overall length 118mm, 4.6 inches.

LBO-400



LBO-401





Kirby Expressor Hook and Loop-Double Ended

LBO-402



5mm wide x 10mm long loop. 10mm wide muscle/expressor hook. Round knurled handle . Overall length 155mm, 6.1 inches.





Paton Double Ended Spatula

LBO-403



6mm wide by 15mm long thin blade with three grooves. 2mm wide by31mm blade with last 7mm angle tip. Flat serrated handle . Overall length 143mm, 5.6 inches.





Nordan Lens Manipulator

LBO-404



Rounded "donut tip" of this instrument minimize chances of posterior capsule rupture. Design also facilities rotation and cracking of the nucleus. Angled tip for increased range of motion when used through a scleral tunnel or clear corneal incision. Round knurled handle. Overall length 120mm, 4.7 inches.

Kellan Nucleus Segmenter

LBO-405



This nucleus segmenter is shaped to fit conveniently into the trough cut into the nucleus by the phacoemulsification tip. With a spreading motion between the phaco tip and the segmenter, the nucleus is cracked. Round knurled handle. Overall length 125mm, 5.0 inches.

Kansas Nucleus Support Spatula

LBO-406



This 3mm wide spatula was designed to slide under the lens nucleus and act as a support for the "Phacofracture" technique using either the kansas Bisection Spatula or the Kansas Trisector. Flat smooth handle. Overall length 132mm, 5.2 inches.



A: Right hand use

B: Left hand use